



Anti-Human IgG(Fcγ fragment specific), AlpSdAbs® VHH

Summary

Code 023-101-001

Immunogen Recombinant Fc region of Human IgG

Host Alpaca pacous

Isotype VHH domain of alpaca IgG2b/2c

Conjugate Unconjugated(6*his tag and one cys were added at the C terminal of the VHH)

Specificity Human IgG(Fcy fragment specific)

Cross-Reactivity Recognizes human IgG Fcy fragment specifically, and reacts with cynomolgus IgG. No Cross-reactivity to rabbit , mouse,

rat, goat IgG

Purity Recombinant Expression and Affinity purified

Concentration 1mg/ml

Formation Liquid, 10mM PBS (pH 7.5), 0.05% sucrose, 0.1% trehalose, 0.01% proclin300 Storage Store at -20 °C(Avoid freeze / thaw cycles), Stable for 12 months at -20 °C

Description

Anti-Human IgG(Fcγ fragment specific), AlpSdAbs® VHH is designed for detecting human IgG Fcγ fragment specifically. Anti-Human IgG(Fcγ fragment specific), AlpSdAbs® VHH is monovalent, recombinant single domain antibodies derived from the variable regions of heavy chain of Alpaca pacous. Based on immunoelectrophoresis and/or ELISA, Anti-Human IgG(Fcγ fragment specific), AlpSdAbs® VHH react with the Fc fragment of human IgG heavy chain but not with the Fab portion of human IgG.

Background

In mammals, antibodies are classified into five main classes or isotypes – IgA, IgD, IgE, IgG and IgM. They are classed according to the heavy chain they contain – alpha, delta, epsilon, gamma or mu respectively. IgG is the most abundant antibody in normal human serum, accounting for 70-85% of the total immunoglobulin pool. Human IgG consists of four human subclasses (IgG1, IgG2, IgG3 and IgG4), and each contains a different heavy chain. The whole IgG molecule possesses both the Fc region and the Fab region, which possessing the epitope-recognition site. The IgG contains two heavy and light chains(kappa or lambda). The heavy chain is about 50 KD and the light chain is about 25 KD. The common IgG is monomeric with a molecular weight of approximately 150 kD.

VHH are single-domain antibodies derived from the variable regions of heavy chain of Camelidae immunoglobulin. The size of VHH is extremely small(<15KDa) compared to other forms of antibody fragment, which significantly increase the permeability of VHH. Thus VHH is considered of great value for research, diagnostics and therapeutics.

Benefits

High lot-to-lot consistency Increased sensitivity and higher affinity Animal-free production

Suggested Working Concentration

ELISA 1:10000-1:50000 WB 1:10000-1:50000

Dilution factors are presented in the form of a range because the optimal dilution is a function of many factors, such as antigen density, permeability, etc. The actual dilution used must be determined empirically.

This product is for research use only and is not approved for use in humans or in clinical

Website: alpvhhs.com E-mail: service@nb-biolab.com Phone: 400-166-9953